# MasterFlux Smart

ETICHETTA:

LABEL:

ÉTIQUETTE:

ETIKETT:

ETIQUETA:

ETYKIETA:

ETIQUETA:

Istruzioni per l'uso	ITALIANO
Instructions for use	ENGLISH
Mode d'emploi	FRANCAIS
Bedienungsanleitung	DEUTSCH
Instrucciones de uso	ESPAÑOL
Instrukcja użycia	POLSKIE
Instruções para o uso	PORTUGUÊS

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---

## CONTENTS

1	INTENDED USE	2
2	NOTES ON SEDATIVE ANALGESIA	3
3	SAFETY	4
3.1	GENERAL SAFETY RULES	4
3.2	IMPORTANT INDICATIONS	5
3.3	SAFETY DEVICES	9
3.4	DEVICE MOVEMENT	9
3.5	- DISPOSAL	10
4	TECHNICAL DATA	10
4.1	AUTONOMY	11
4.2	ELECTRONIC AUTONOMY MANAGEMENT (OPTIONAL)	12
5	COMPONENTS	13
5.1	COMPONENTS SUPPLIED	13
5.2	CONTROL PANEL	14
5.3	TOUCH SCREEN DISPLAY	15
5.3.1	MAIN MENU	15
5.3.2	SEDATION MENU	18
5.4	BATTERY	19
5.5	CYLINDERS	20
6	SEDATION	21
6.1	CIRCUIT AND ANAESTHESIA BAG CONNECTION	21
6.2	CIRCUITS AND MASKS	22
6.3	INDICATIONS FOR CONNECTING TO THE MAINS, SWITCHING ON THE DEVICE AND USB CONNECTION	23
6.4	PASSIVE AND ACTIVE SUCTION CONNECTION	24
6.5	PRECAUTIONARY MEASURES TO PREVENT EXPOSURE TO N2O	26
6.6	START AND END OF SEDATION	27
6.7	SEDATION RECORDING, ARCHIVING AND EXPORTING	29
7	ALARMS	30
7.1	WARNINGS	30
7.2	- ERRORS	31
7.3	WI-FI COMMUNICATION	31
8	MAINTENANCE	32
8.1	ROUTINE MAINTENANCE - CLEANING	32
8.2	ROUTINE MAINTENANCE - CYLINDER REPLACEMENT	32
8.3	EXTRAORDINARY MAINTENANCE	33
8.4	PROBLEMS - FAILURE TO SEDATE THE PATIENT	33
9	SUPPORT	34

## 1 INTENDED USE

MASTER FLUX SMART is a device for mixing and administering medical O<sub>2</sub> and N<sub>2</sub>O to the patient, used in conscious sedation / sedative analgesia techniques.

### USE

The concentration values in the administered mixture are:

O<sub>2</sub> ≥ 30%

N<sub>2</sub>O ≤ 70 %

The administration is performed with a nasal mask, the patient must inhale the mixture by actively breathing through the nose.

Conscious sedation is indicated for all patients in the dental practice, both adults and children.

Conscious sedation allows patients to undergo dental treatment while limiting the fear, anxiety and sense of discomfort they may experience during dental treatment.

### SAFETY MARKINGS AND SYMBOLS USED

	CE marking proving the conformity of the Product with Directive 93/42EEC and subsequent amendments and supplements
	Type B applied part - is part of an electromedical device that in normal use necessarily comes into physical contact with the patient. The leakage current on this part must be lower than the limits set by the standard.
	Date of manufacture (month and year)
	Fuses used in the device
	It is mandatory to read the manual
	Disposal
	Model
	Serial number
	Mixed gas outlet + exhaust vac
	Passive exhaust + active exhaust
	Warning sticker + ON/OFF

The Declaration of Conformity is attached to the device

## 2 NOTES ON SEDATIVE ANALGESIA

The gas sources are nitrous oxide (N<sub>2</sub>O) and oxygen (O<sub>2</sub>), both for medical use.

The gas mixture delivered enters the administration circuit, is inhaled by the patient, and then discharged outside the practice.

There are two one-way valves on the mask to ensure proper gas flow.

### PRINCIPLE OF ANALGESIA

Nitrous oxide is a laughing gas, it is an inert, colourless, sweet-tasting, chemically stable, non-flammable and explosive gas, but it can fuel combustion.

It is a non-toxic narcotic for inhalation, which does not stimulate the respiratory tract and does not combine with haemoglobin.

It takes effect quickly after being inhaled into the body. It begins its analgesic activity within 30-45 seconds of inhalation.

Normal breathing allows 90% of nitrous oxide to be eliminated, 10% evaporates from the patient's skin more slowly

**SEDATION SHOULD ALWAYS BE MONITORED USING A PULSE OXIMETER WITH CUSTOMISED ALARM THRESHOLDS.**

Vital parameters, pulse and oxygen saturation show whether the patient is in good condition and whether sedation has the desired effect.

Normally, during sedation, oxygen saturation increases, the pulse slows down and the patient relaxes.

When a patient is sedated with e.g. 50% nitrous oxide, the oxygen flow is 50%, the patient breathes more oxygen than normal (oxygen in the air ~ 21%).

**If the oxygen saturation falls below the set value, the pulse oximeter emits an alarm.**

If oxygen saturation falls below 85%, for example, nitrous oxide sedation must be discontinued immediately, as there is a risk that nitrous oxide and oxygen have been exchanged and the mixture administered contains 80% nitrous oxide.

**If beats rise or fall above the set value, the pulse oximeter will sound an alarm.**

It could indicate an overdose or too deep sedation. In this case, reduce or set to 0 the flow of nitrous oxide and administer pure oxygen to the patient.

If the saturation indicated by the pulse oximeter falls during the administration of pure oxygen, switch off the unit immediately and contact service.

## 3 SAFETY

### 3.1 GENERAL SAFETY RULES

	<b>The MASTER FLUX SMART must only be used by specially trained medical personnel. (*)</b>
	<b>Use is only permitted in the areas indicated in the intended use</b>

The use of nitrous oxide (N<sub>2</sub>O) in the out-of-hospital setting should not be combined with the use of other anaesthetic drugs, hypnotics, sedatives or major painkillers except in the presence of a resuscitating anaesthetist. On the other hand, the use of nitrous oxide may precede or follow the infiltration of local anaesthetics, such as in dentistry.

	<b>Please read this user manual carefully. The user manual is an integral part of the equipment. Keep it close to the unit for reference at all times.</b>
---	--

**(\*) The text is valid for devices sold in Italy.**

TECNO-GAZ recommends the application of this text for all users of the device, integrating it with the regulations in force in the country of use.

AIFA - DDL guidance on the extra-hospital use of nitrous oxide

The use of nitrous oxide (N<sub>2</sub>O) is only permitted through specific equipment that prevents the delivery of hypoxic mixtures and only in specific environments with adequate air circulation.

The drug is to be used in an out-of-hospital setting at a maximum concentration of 50% oxygen (i.e. no more than 50% nitrous oxide and at least 50% oxygen).

The only exception is conscious sedation practised by experienced sedationists who, with the patient definitely conscious (as documented by the opening of the mouth), can set up, for short periods of time, a delivery at concentrations of up to 70%.

Nitrous oxide has been widely used for decades in non-hospital settings, particularly in dentistry, as an anxiolytic analgesic. No serious complications secondary to its use as an anxiolytic analgesic in out-of-hospital settings have been described to date.

Nitrous oxide must only be used in the presence of medical or dental personnel who are trained in cardiopulmonary resuscitation (BLS-D) and who are aware that A) in the event of loss of consciousness, nitrous oxide should be immediately discontinued and B) in the event of apnoea, the patient's airway should be immediately opened (e.g. with the threefold manoeuvre of extending the head, subluxating the jaw and opening the mouth) and, if necessary, ventilation support should be provided (e.g. with a self-expanding balloon and a face mask).

## 3.2 IMPORTANT INDICATIONS

### PATIENT SELECTION

Conscious sedation with nitrous oxide is recommended for all patients who visit the dentist with discomfort or fear, as it allows them to experience the treatment in a relaxed manner.

- Scared, frightened patients
- Patients with a strong gag reflex
- Patients in whom insufficient local anaesthesia can be achieved
- Cooperative children in need of long and extensive treatment
- Children with limited cooperation but less need for treatment

### CONSCIOUS SEDATION IS NOT INDICATED FOR PATIENTS:

- IN PREGNANCY

(even dental office employees cannot be exposed to nitrous oxide)

- IN DRUG-ADDICTED PATIENTS
- IN THE PRESENCE OF SEVERE PULMONARY INFECTIONS
- SEVERE MENTAL OR BEHAVIOURAL AND PERSONALITY DISORDERS
- PROBLEMS OR LIMITATIONS IN NASAL BREATHING (RHINITIS, SINUSITIS)
- REDUCED COMMUNICATION SKILLS
- DISTURBANCE OF MIDDLE EAR VENTILATION
- SERIOUS SYSTEMIC DISEASES (SM, ILEUM)
- HISTORY OF SURGERY WITH INTRAOCULAR GAS INJECTION
- PATIENTS WITH ANESTHESIOLOGICAL EVALUATION EQUAL TO OR ABOVE 3 (See table below)
- VITAMIN B12 AND FOLIC ACID DEFICIENCY OR RELATED DISORDERS

Class ASA	Definition	Suitability for nitrous oxide sedation
1	Normal healthy patients	Excellent
2	Patients with mild systemic diseases	In general, good
3	Patients with severe systemic diseases	Do not sedate with laughing gas in the practice
4	Patients with severe systemic diseases in constant danger of dying	Not good risks prevail
5	Dying patients, who would probably not survive without intervention	Not suitable
6	Brain-dead patients to be harvested for organ donation	Not suitable

ABSOLUTE CONTRAINDICATIONS	
Ileo Pneumothorax Pneumopericardium Pneumomediastinum Enterothorax Lung cyst Middle ear surgery Middle ear inflammation Gas embolism Cranio-cerebral trauma Intraocular air bubble Perforating eye injuries Ocular fund operations	As nitrous oxide penetrates air-filled cavities and spaces depending on the concentration gradient, it causes a local increase in volume
Psychotropic drugs Oral sedatives Drug addiction Alcoholism	Interactions with the substances listed here are not predictable.
Vitamin-B12/folic acid deficiency or associated disorders	Nitrous oxide oxidises the cobalt in vitamin B12, interfering with erythropoiesis.
Limited nasal respiration	Prerequisite for effective sedation continuous nasal breathing
COPD chronic obstructive pulmonary disease	Since in patients with COPD (chronic obstructive pulmonary disease) the respiratory stimulus driven by the blood oxygen concentration, the administration of oxygen during the sedation could suppress this stimulus.
Increased intracranial pressure Alteration of the blood-brain barrier	In these patients, extravasation of fluids can lead to increased intracranial pressure.
Left ventricular failure Right heart load	Nitrous oxide has an effect on the cardiodepressive (Follow ASA classification).
Bleomycin therapy	Up to 18% of patients treated with bleomycin develop pulmonary fibrosis.
Pregnancy	Embryonic damage cannot be excluded.
Personality disorders/severe behavioural abnormalities psychosis/phobias	The action of nitrous oxide may lead to a worsening of symptoms in these patients.
Cystic fibrosis (mucoviscidosis)	Mucoviscidosis involves chronic inflammation of the airways.
Drug addiction or withdrawal	Nitrous oxide can cause undesirable behaviour.
Claustrophobia	Patients often do not tolerate the nasal mask.
RELATIVE CONTRAINDICATIONS	
Disabled patients	Consultation with the treating physician is necessary as pharmacotherapy is often complex. It must be ensured that the patient breathes through the nose.
Seniors	Even in geriatric patients it is imperative to consult with the attending physician.

(Taken from "Conscious sedation with nitrous oxide in the dental surgery - Wolfgang Lüder)

Nitrous oxide diffuses into the cavities, so it should not be inhaled after a middle ear infection, intestinal occlusion, ophthalmic surgery with intraocular gas injection, etc.

**Always take an anamnesis and follow what was learned during the training course.**

## OVERDOSE

Excessive dosage of nitrous oxide may cause dizziness, acceleration of pulse and malaise. Laughing, tiredness, delayed reaction to questions, difficulty in keeping the mouth open, feeling of "flying", sweating and wrinkled forehead to the point of fear, migraine, nausea, vomiting, nightmares in case of significant overdose.

The operator must learn how to titrate nitrous oxide for the individual patient.

In the event of an overdose, cut off the nitrous oxide supply or press the FLUSH button.

## GENERAL WARNINGS

	<b>The law restricts the sale of this device to MEDICAL AND DENTAL PRACTITIONERS only.</b>
	<b>Use only MEDICAL Oxygen and Nitrous Oxide cylinders</b>
	<b>Warning: The device is not designed for home use.</b>
	<b>The device must not be used in conjunction with anaesthesia equipment or as a system for the administration of anaesthesia.</b>
	<b>The device only administers medication for the conscious sedation technique.</b>
	<b>Before each use, check the correct connection of the supply gases.</b>
	<b>Contraindications to the use of the medicine are given in the safety data sheet enclosed with the medicine. If not present, ask your gas supplier.</b>
	<b>Slowly open the valves of the oxygen and nitrous oxide cylinders.</b>
	<b>Only use accessories supplied by TECNO-GAZ to ensure the safety of the patient and the operator.</b>
	<b>Modifications to the device are not permitted</b>
	<b>Unauthorised repair, use of non-original spare parts, alteration or abuse of this device will invalidate the warranty and in case of damage to things or persons Tecno-gaz will not be held responsible.</b>

	<b>Sedation should always be monitored using a pulse oximeter with customised alarm thresholds. Always keep spare batteries available for the pulse oximeter.</b>	
	<b>Accessory not supplied (Art. 11260)</b>	
	<b>O2 saturation levels</b>	<b>Degree of hypoxia</b>
	<b>≥ 95%</b>	<b>No hypoxia</b>
	<b>90-94%</b>	<b>Mild hypoxia</b>
<b>85-89%</b>	<b>Moderate hypoxia</b>	
<b>84 %</b>	<b>Severe hypoxia</b>	

**ELECTRICAL WARNINGS**

	<b>Connect the device to a socket with protective earth</b>
	<b>Noise in the surrounding environment can generate interference on the device</b>
	<b>Interference between the device and other devices in the same working environment is possible</b>

### 3.3 SAFETY DEVICES

The device is equipped with safety, pneumatic and mechanical systems:

If the oxygen (O<sub>2</sub>) pressure is reduced, the nitrous oxide (N<sub>2</sub>O) supply is automatically interrupted.

A built-in mixing system regulates the N<sub>2</sub>O in such a way that a minimum of 30% O<sub>2</sub> is always present in the gas mixture.

If the flow of the gas mixture is interrupted, the patient sucks in air from outside via a safety valve.

Should the flow change direction back towards the device, a one-way valve protects the patient from breathing in exhaled gases and accumulating carbon dioxide (CO<sub>2</sub>).

The FLUSH button, when pressed, interrupts the flow of N<sub>2</sub>O and generates a high flow of pure oxygen.

THE DEVICE IS EQUIPPED WITH ELECTRONIC SAFETY DEVICES THAT MONITOR AND SIGNAL THE ANOMALIES IN THE SYSTEM.

	<p>IT IS ESSENTIAL THAT THE OPERATOR DURING ADMINISTERING:</p> <ul style="list-style-type: none"> <li>• IS ALWAYS PRESENT AND KEEPS AN EYE ON THE PATIENT'S VITAL FUNCTIONS</li> <li>• KEEPS FLOWS UNDER CONTROL</li> <li>• CHECKS THE BAG FOR OPTIMAL PATIENT VENTILATION</li> </ul>
	<p>EXCESSIVE BRIGHTNESS MAY MAKE IT DIFFICULT TO READ THE DISPLAY. CHECK THAT THE % AND FLOW DELIVERED ARE AS DESIRED.</p>

### 3.4 DEVICE MOVEMENT

To move the device

- Check the correct positioning of accessories (circuits, masks, etc.)
- Check that the power cable is not connected and positioned correctly.
- Check that the cylinders are closed.


**Unlock the 4 wheels and use only the handles on the sides of the cabinet for movement.**

 <p>Do not move the cabinet without using the handles provided.</p>	 <p>Do not climb on top of the cabinet.</p>	 <p>Do not sit or place weights on the cabinet.</p>
--	--	--

**Once the movement is complete, secure the device by locking the wheels.**

### 3.5 DISPOSAL

Disposal of the packaging, device and individual components must be in accordance with the regulations in force in the country of use for the treatment and disposal of potentially bio-contaminated waste. Do not dispose of packaging material in the environment.

The Declaration of Conformity is attached to the device

## 4 TECHNICAL DATA

<b>CLASSIFICATION</b>	
<i>Classification of the medical device according to Directive 93/42/EEC updated with the changes made by the new Directive 2007/47/EC (Annex IX)</i>	<i>IIb (in application of rule 11) Short-term duration Non-invasive medical device. Active medical device</i>
<b>SAFETY CATEGORY</b>	
<i>Safety class</i>	<i>Type I</i>
<i>Level of protection against electric shock</i>	<i>Applied part type B</i>
<b>O<sub>2</sub> AND N<sub>2</sub>O SPECIFICATIONS</b>	
<i>Supply gas flow rate (l/min)</i>	<i>~ 100 LPM</i>
<i>Supply gas pressure (bar / PSI)</i>	<i>3.5 ÷ 6 bar - 50 ÷ 87 PSI</i>
<i>Maximum device flow rate (l/min)</i>	<i>12 NI/min</i>
<i>Maximum gas supply pressure (bar / PSI)</i>	<i>6 bar - 87 PSI</i>
<b>ELECTRICAL SPECIFICATIONS</b>	
<i>External mains power supply (Power supply input)</i>	<i>AC 100 ÷ 240 V +/-10% 50 ÷ 60 Hz</i>
<i>Built-in battery: Pb-Ca battery</i>	<i>Capacity: 7Ah; Rated voltage 12VDC</i>
<i>Power consumption</i>	<i>40VA</i>
<i>Battery fuse 12VAC</i>	<i>F2A (5 x 20 mm)</i>
<i>Mains fuses</i>	<i>F1A (5 x 20 mm)</i>

<b>ENVIRONMENTAL SPECIFICATIONS: WORKING CONDITIONS</b>	
<i>Operating temperature</i>	5 ÷ 35 °C
<i>Relative humidity (RH)</i>	≤90% (non-condensing)
<i>Atmospheric pressure</i>	76 ÷ 106 kPa
<b>ENVIRONMENTAL SPECIFICATIONS: TRANSPORT AND STORAGE CONDITIONS</b>	
<i>Transport and storage temperature</i>	-5 ÷ 40 °C
<i>Relative humidity (RH)</i>	≤90% (non-condensing)
<i>Atmospheric pressure</i>	76 ÷ 106 kPa
<b>WEIGHT AND DIMENSIONS:</b>	
<i>Weight without cylinders</i>	47 kg
<i>Overall dimensions of the device H; W; D</i>	H1230 x D590 x W410 [mm]
<b>CYLINDERS:</b>	
<i>Maximum height</i>	950 mm
<i>Maximum diameter</i>	140 mm


#### 4.1 AUTONOMY

<b>EXAMPLE OF CYLINDER AUTONOMY</b>	
<i>MEDICAL OXYGEN BOTTLE WITH CAPACITY 10 l. LOADED AT A PRESSURE OF 200 BAR</i>	
<i>AVAILABLE VOLUME</i>	2000 NI
<i>NITROUS OXIDE CYLINDER WITH 10 l CAPACITY. LOADED WITH 7 kg OF LIQUID GAS</i>	
<i>AVAILABLE VOLUME</i>	4600 NI
<b>SETTING:</b>	
<i>DELIVERY</i>	10 NI/min.
<i>MIXING</i>	50 %
<b>AUTONOMY:</b>	
<i>OXYGEN</i>	400 min.
<i>NITROUS OXIDE</i>	950 min.
	Based on average gas consumption TECNO-GAZ recommends replacing the nitrous oxide cylinder every 2 Oxygen cylinders.

## 4.2 ELECTRONIC AUTONOMY MANAGEMENT (OPTIONAL)

It is also possible to manage the autonomy automatically, but it is necessary to reset the cylinders each time a full cylinder is fitted.

From the MENU press VIEW, set the cylinder volume and press on the DATE to reset.

The device automatically estimates the autonomy and displays the REMAINING MINUTES (5) in delivery and the NORMAL LITRES available if there is no delivery, the estimate is made by subtracting the litres delivered from the volume of gas available.

i	<b>When the estimated autonomy of a gas is reduced to less than 30 minutes, the autonomy box will turn yellow and flash</b>
i	<b>Tecno-gaz does not accept responsibility for incorrect assessment of autonomy.</b>

## 5 COMPONENTS

The device consists of a cabinet with 4 swivel wheels with brakes.

At the top of the front is the control panel, which is tilted for easy reading.

Underneath the control panel there is the delivery block, where the tubes of the gas delivery circuit are connected to the patient and where the anaesthesia bag is connected, this space can be used for storing the accessories.

At the rear there is a storage area designed to accommodate two cylinders with a closing door.

On the right side, there is a handle at the top for moving the device.

At the bottom there are connectors for ACTIVE and PASSIVE discharge of EXHALED GASES.

On the left side, there is a handle for moving the device and a USB socket.

At the bottom, there is a main switch and a plug for connecting the device to the power supply.

Inside the device there is a battery that allows the use of the device even without connection to the mains.

### 5.1 COMPONENTS SUPPLIED

ITEM	ACCESSORIES KIT	N <sub>2</sub> O REDUCER	O <sub>2</sub> REDUCER	CIRCUITS AND MASKS SUPPLIED
<b>MF001IMF</b> Master Flux Smart with Reducers with UNI Oxygen and Nitrogen Oxide connections	2IMFA0030	MF011ZMF	MF847ZMF	MF820ZMF Administration circuit
<b>MF002IMF</b> Master Flux Smart with Reducers with AFNOR Oxygen and Nitrogen Oxide connections	2IMFA0051	MK005ZMF	MK012ZMF	MF800IMF Active discharge circuit SMFA197 Passive discharge circuit
<b>MF003IMF</b> Master Flux Smart with Reducers with DIN Oxygen and Nitrogen Oxide connections	2IMFA0052	MF011ZMF	MF845ZMF	MF835ZMF Intelliflux 1 MF836ZMF Intelliflux 2 MF837ZMF Intelliflux 3
<b>MF004IMF</b> Master Flux Smart without reducers kit	2IMFA0053	-	-	

## 5.2 CONTROL PANEL

The control panel consists of

### 1 - DELIVERY KNOB

Turning clockwise increases the flow rate, and allows the flow rate to be set.

### 2 - DOSING KNOB

Turning clockwise increases the % of nitrous oxide delivered (from 0 to 70%).

### 3 - FLUSH BUTTON

Pressing stops the nitrous oxide supply and provides a high flow of oxygen.


### 4 - TOUCH SCREEN DISPLAY

Used to manage settings, store sedation and display parameters.

### 5 - BAR LED

The BAR LED indicates the status of the machine according to the colour:

GREEN	device switched on and connected to the mains
BLUE	device in delivery phase
YELLOW	device with WARNING signalled
RED	device in ERROR
OFF	device powered by the internal battery or switched off

### 5.3 TOUCH SCREEN DISPLAY

#### 5.3.1 MAIN MENU

At the top of the display there is some general information that is offered in all screens:





**DOCTOR (rosette symbol indicates ADMINISTRATOR)**





DATE, TIME, WIFI (if present), % BATTERY CHARGE.

When the device is started for the first time, the administrator is registered, who can enter other doctors. The doctor who uses the device can then log in directly by entering the password set during registration. At this point there is the possibility to choose one of the three available menus:

<p><b>DOCTOR</b> - Only the administrator can give access to new doctors. While administrator and doctors can change their personal password and remove their access credentials.</p>	
<p><b>PATIENT</b> - It is used to enter a new patient by filling out the appropriate form. Select a patient or perform anonymous sedation. Choose whether to record sedation. From here the SEDATION MENU (5.3.2) can be displayed</p>	
<p><b>MENU</b> - consisting of the following items:</p>	

<p><b>DOCUMENTS</b> allows to:</p> <ul style="list-style-type: none"> <li>- display the last sedations stored (max. 99)</li> <li>- export to USB stick</li> </ul>	
<p><b>DISPLAY</b> is used to:</p> <ul style="list-style-type: none"> <li>- set the volume of cylinders used and do a cylinder reset</li> <li>- set the autonomy display (see 4.1)</li> <li>- set the suction flow display (in case of active suction)</li> </ul>	
<p><b>CONNECTION</b></p> <p>Allows you to set up the WIFI connection. You need to enter the network name and password.</p> <p>WI-FI network settings:</p> <ul style="list-style-type: none"> <li>- presence of the password,</li> <li>- 2.4 GHz band,</li> <li>- internet protocol IPV4,</li> <li>- WPA and WPA2 security protocol</li> </ul> <p><b>The device does not support 5 GHz band communication.</b></p>	
<p><b>SERVICE</b></p> <ul style="list-style-type: none"> <li>- INFORMATION displays information about the device (e.g. serial number)</li> <li>- PERIPHERAL CONTROL is used to control and view individual components (access for authorised technicians only)</li> <li>- SERVICING AND MAINTENANCE is used to record (for authorised technicians only) and display the maintenance carried out.</li> </ul>	

<p><b>PATIENT MANAGEMENT</b> is used to select patients, edit patient master data, delete and create new patients.</p>	
<p><b>VOLUME</b> is used to set the intensity of the sound signal emitted when the touch screen is pressed.</p>	
<p><b>LANGUAGE</b> is used to set the language.</p>	
<p><b>DATE AND TIME</b> is used to set the date and time and the time zone.</p>	

## 5.3.2 SEDATION MENU

Screen description:

1 - PATIENT

2 - TOTAL FLOW (by pressing on the total flow it is possible to display the individual flows of O<sub>2</sub> and N<sub>2</sub>O in position 4)

3 - N<sub>2</sub>O / O<sub>2</sub> PERCENTAGE (pressing it displays the O<sub>2</sub> and N<sub>2</sub>O % alternately)

4 - O<sub>2</sub> / N<sub>2</sub>O FLOW

5 - AUTONOMY, the available volume in NORMAL LITRES is displayed. If the device is delivering, the remaining time in MINUTES is displayed. It is possible to choose whether to display the autonomy, using the appropriate command (press menu → view and select the appropriate item).

6 - VACUUM, the SUCTION FLOW is displayed. It is possible to choose whether to display the suction flow by using the appropriate command (press menu → display and select the appropriate item). The flow is only displayed when using the active suction circuit.

7 - BASELINE the baseline values (total flow and N<sub>2</sub>O %) are displayed, if previously set. To set them press the Save Baseline icon during sedation, this saves the flow of O<sub>2</sub> and N<sub>2</sub>O that the device is delivering. It is also possible to set/change it directly on the patient form. It is not possible to set O<sub>2</sub> values lower than 30%.

8 - % N<sub>2</sub>O / % O<sub>2</sub> BAR is a direct display of the percentage of gas delivered.

9 - DISPLAY / RECORDING (red dot when recording).

Recording starts when more than 0.5 nl/min starts to be delivered and ends when gas delivery drops below 0.5 nl/min.

## 5.4 BATTERY

The device is equipped with a BATTERY that guarantees operation even in the event of a mains power failure.

The charge is displayed on the screen

The representation is:

WHITE from 100% to 21%

YELLOW from 20% to 11%

RED from 10% to 0%

The autonomy is approximately 6 hours operating time (on delivery), but varies depending on factors such as % of charge and battery life (charge/discharge cycles).

**To charge the battery, the device must be switched on.**

The recommended time for a full charge is 24 hours.

When the battery charge is low, the device indicates on the display, and with an acoustic message, to switch off the device or connect it to the power supply, at this stage the device cannot be used.

### **Warning:**

When charging the battery the charge indicator rises progressively, if the device is restarted during charging, the indication may overestimate the actual battery charge. In this case, at the start of the cycle, the charge level may decrease rapidly until the actual charge value is reached.

**SWITCH OFF OR CONNECT THE DEVICE TO THE MAINS WHEN IT IS DISCHARGED TO AVOID DAMAGING THE BATTERY**

**WE RECOMMEND REPLACING THE BATTERY EVERY 18 MONTHS.  
THE INTERVENTION MUST BE CARRIED OUT BY AUTHORIZED TECHNICIANS.**

## 5.5 CYLINDERS

The cylinders are located in the rear compartment. 5 - 10 cylinders can be fitted, with the limitation that the maximum diameter of the cylinder must not exceed 140 mm and the height must not exceed 900 mm; they must be fixed to the wall of the cabinet by means of the special strap provided.

The cylinders must be connected to the device via the pressure reducer provided.

The connection is made without the use of tools, the cylinder-reducer connections are gas-specific while the reducer-device connection for O<sub>2</sub> has a right-hand thread and for N<sub>2</sub>O a left-hand thread.

The reducers may cool down and become covered with frost, but this does not prevent operation.

**Carefully follow the instructions in the instruction manual of the pressure reducer.**

Main notes for the use of reducers:

- A. Check that hands are clean before each intervention**
- B. Check the integrity of the reducer**
- C. Check the integrity of the gasket between the reducer and the cylinder**
- D. Do not use oils, grease, solvents or detergents of any kind for cleaning**
- E. Open valves slowly**

Connect the outlet of the reducer via the appropriate fittings to the hoses in the cylinder compartment.

**Pay attention to the reducer-device connection:**

- OXYGEN has a RIGHT-hand thread**
- NITROGEN OXIDE has a LEFT-hand thread**

**At the end of the day, store the cylinders in a ventilated storage area in accordance with the laws in force**


## 6 SEDATION

### 6.1 CIRCUIT AND ANAESTHESIA BAG CONNECTION

Below the control panel there is the delivery module, which consists of three 22 mm diameter male connectors.


- A. Two connectors on the left side: a vertical one dedicated to the connection of the anaesthesia bag and a 45° one for the connection of the gas delivery tube to the patient.
- B. On the right side there is a single connector dedicated to the connection of the gas evacuation tube from the patient.

## 6.2 CIRCUITS AND MASKS

**Follow the instructions dedicated to the circuits in the use and maintenance manual.**

### MASK

①	Support zone for the patient's forehead
②	Gas inlet with non-return valve
③	Gas outlet with non-return valve + arrow printed on PC to indicate flow direction
④	Velcro elastic band with plastic hooks

### MASK SIZES

MF835ZMF	MASK 1	0-4 years
MF836ZMF	MASK 2	4-10 years
MF837ZMF	MASK 3	>10 years
MF838ZMF	MASK 4	

### CIRCUIT CONSISTING OF

NO. 2 SUPPLY AND DISCHARGE PIPES

PIPE SUPPORT CLIP AND VELCRO BAND

### 6.3 INDICATIONS FOR CONNECTING TO THE MAINS, SWITCHING ON THE DEVICE AND USB CONNECTION

1 - The mains cable must be connected on the right-hand side of the device. The socket (1) has two protective fuses.

2 - Above the socket there is the switch (2) which starts the device.

3 - At the top near the handle there is the USB socket - the USB socket allows to download the reports of the sedations performed and if necessary to update the firmware of the device.

The USB socket when not in use:

**Must be protected by the plastic cover provided.**


**Beware of electrostatic discharges that the operator may generate.**

**The single-pole switch does not allow the device to be disconnected from the mains.**

**TO DISCONNECT THE DEVICE FROM THE MAINS, THE POWER CABLE MUST BE UNPLUGGED.**

## 6.4 PASSIVE AND ACTIVE SUCTION CONNECTION

To reduce the exposure of operators to nitrous oxide, it is necessary to take the exhaled gases outside the practice.

Connect the hose coming from the patient's mask to the appropriate connection B, par 6.1.

Choose the type of evacuation and proceed as follows:

### ACTIVE SUCTION

Connect the pipe (included in the supply) to the appropriate fitting pos 3.


Adjust the aluminium valve with adjustment ring nut to limit the suction, open the valve to reduce the suction, close the valve to increase the suction (fine adjustment – pos. 4)

Adjust the suction terminal to close and open the suction (coarse adjustment), connect this terminal to a vacuum source.

#### Attention the vacuum must:

- discharge outside the practice
- be suitable for anaesthetic gases suction

The flow of suctioned gas is visible on the display (see 5.3 SEDATION MENU).

**BY VARYING THE FLOW DELIVERED, THE SUCTION FLOW MUST BE ADJUSTED.**


#### Warning:

**THE TOTAL DELIVERED FLOW (5) MUST BE ALIGNED WITH THE SUCTION FLOW (6).**


### PASSIVE SUCTION

Unscrew the red cap in position 7 and, in its place, screw the supplied hose connector. Place the other end, where the filter is positioned, outside the cabinet in a ventilated environment.

**ACTIVE SUCTION WITHOUT SUCTION FLOW MEASUREMENT**


Remove hose connector A  
Remove hose connector B and fit hose connector (A) removed earlier

Remove the red cap.  
Screw the modified hose connection to the hose fitting position 7.

Follow the instructions above:  
"SUCTION ACTIVE".

We also recommend going into the MENU, enter VIEW and uncheck:  
"set suction flow display".  
See par. 5.3.1

## 6.5 PRECAUTIONARY MEASURES TO PREVENT EXPOSURE TO N<sub>2</sub>O

Environmental concentration limits of N<sub>2</sub>O. Maximum environmental concentrations of N<sub>2</sub>O must be measured in parts per million (ppm).

The value must not exceed the imposed limits.

Compliance with the standards and precautionary measures listed below can reduce N<sub>2</sub>O concentrations to values well below the required limits.

- the morning before use, check the rubber hoses, connections, any other tubing, the balloon and their tightness and replace them if damaged or at the proposed replacement date;
- have masks of varying sizes so that they fit as comfortably and securely as possible on the patient's nose and face. Leaks through masks are the most important sources of N<sub>2</sub>O in the area where the dentist and assistant work.
- make sure that the exhaled gas suction system is set correctly. Active suction is the most effective method of reducing environmental N<sub>2</sub>O concentration.
- The patient should be asked to restrict speech and not to breathe through the mouth after the mask has been applied to the nose.
- in children, crying, patient movements, the procedures used including the application of a topical anaesthetic and the performance of loco-regional anaesthesia and the application of a mouth opener, when the patient talks, laughs, cries or is uncooperative may be additional factors for pollution.
- if available, activate an environmental ventilation system; staff should not stand between the patient and the ducting that directs ambient air to the ventilation system outlet. Doors and windows should be open.

At the end of the procedure the patient should be asked to breathe 100% O<sub>2</sub> through the mask for at least five minutes. This time interval is necessary to eliminate N<sub>2</sub>O residues.

In paediatric dentistry, reluctant patients who do not tolerate the nasal mask and patients with medical conditions that contraindicate the use of N<sub>2</sub>O (obstructive airway disease, emotional disorders, drug dependence) should be treated with other techniques.

## 6.6 START AND END OF SEDATION

Check the integrity of the components before each use

Before starting treatment, carry out a clinical, psychological and physical assessment of the patient (blood pressure, radial pulse and respiration).

The device does not monitor the patient's physiological parameters, it is necessary to monitor the patient using a pulse oximeter with customised alarm thresholds.

Place the cabinet in a flat, stable, comfortable position visible to the operator. Lock the wheels.

Allow proper air exchange by simply keeping a window open.

Slowly open the valve of the oxygen cylinder, check the value on the pressure gauge, do the same for nitrous oxide.

Check that the anaesthesia bag is correctly installed, connect the delivery circuit and the exhaled gas circuit according to the instructions in the manual enclosed with the components.

Switch on the device and enter the doctor's credentials.

Enter the patient's details if not already present or proceed with anonymous sedation.

Confirm, if desired, the sedation registration.

**At this stage do a functional check, close the oxygen cylinder and check that the device signals a lack and that the flow of nitrous oxide is blocked, if the operation is correct reopen the cylinder and proceed with sedation, if there are anomalies stop sedation close the cylinders and contact technical support.**

Choose the most suitable mask for the patient and apply it to the patient, taking care that the edge has a good fit on the patient's face.

Start delivering by turning the knob (O<sub>2</sub>) on the left-hand side of the control panel.

The device starts recording.

Tell the patient to breathe through the nose only.

Turn the knob until the optimal O<sub>2</sub> flow for the patient is achieved, the correct flow is checked by looking at the anaesthesia bag, it should remain open but not inflated **"half full / half empty"** and should follow the patient's breathing.

THERAPY SPECIFICATIONS FOR CHILDREN:	THERAPY SPECIFICATIONS FOR ADULTS:
Set a flow rate of 4/5 l/min and allow to breathe in for ~1 minute.	Set a flow rate of 6-8 l/min and allow to breathe in for ~1 minute.

If passive discharge is used, there are no special warnings.

If active suction is used, it is important to adjust the suction flow effectively (see 6.3).

It is important to adjust the active suction correctly, if the suction flow is high the patient does not enter sedation and the anaesthesia bag tends to close.  
If the suction is too low, the anaesthesia bag tends to swell.  
In both cases, the patient's breathing is not followed.

Once the correct flow has been determined (depending on the patient's pulmonary capacity) it is necessary to progressively increase the % of N<sub>2</sub>O starting from 10% and increasing in steps of 5% using the right hand knob, wait at least one minute before increasing the dosage, until the patient is sedated.

**The values, flow and %N<sub>2</sub>O are characteristic of the patient and are called BASELINE or TITRATION. It is possible to store this data in the registry by pressing the appropriate button below the patient's name.**

The device will display these values every time we display the patient in the sedation menu and in the patient registry.

It will be possible to modify the saved data by entering the patient data.

**Obtaining the correct baseline causes a state of relaxation and well-being in the patient.**

**The baseline ensures the minimum amount of drug (N<sub>2</sub>O) required for the intervention, and also allows rapid removal of the drug from the body.**

The intervention performed on the patient (sedated but awake) will be accompanied and remain in the memory with the perception of a positive experience.

**BEFORE TERMINATING THE DELIVERY, TURN THE %N<sub>2</sub>O KNOB TO THE 0 POSITION AND GIVE THE PATIENT OXYGEN ONLY FOR 5 MINUTES.**

The device will terminate recording when the flow is reset.

	<p>If the patient experiences symptoms of malaise (nausea, dizziness, headache), STOP THE THERAPY IMMEDIATELY. There is a "FLUSH" button on the device. Pressing this button will fill the balloon with oxygen. By squeezing the balloon we administer high quantities of pure oxygen to the patient.</p>
	<p>CONSTANTLY MONITOR THE FLOWS SHOWN ON THE DISPLAY DURING THERAPY</p>
	<p>BEFORE LEAVING THE PATIENT, EVEN FOR SHORT PERIODS, STOP SEDATION.</p>
	<p>THE PERCENTAGE OF NITROUS OXIDE ADMINISTERED TO THE PATIENT SHOULD BE AS LOW AS POSSIBLE TO ACHIEVE GOOD SEDATION.</p>
	<p>Before discharging the patient, keep him/her under observation (waiting room) for about 5-10 minutes. The doctor may deem it necessary to increase the observation time.</p>
	<p>After using the device, it is IMPORTANT to close the cylinder valves.</p>

## 6.7 SEDATION RECORDING, ARCHIVING AND EXPORTING

The device records sedation, if the operator chooses this option.

Recording starts when medical gases are delivered and ends when delivering stops.

Warnings / errors that may occur during sedation are recorded

**The following data is recorded:**

- Doctor
- Patient
- Date and time
- Device
- Device serial number

**In relation to time the following is recorded:**

- Flow
- % of nitrous oxide
- Warnings / Errors

**Recording takes place every time a parameter (flow, % N<sub>2</sub>O) is changed.**

**If sedation is interrupted and resumed after a few seconds there will be two recorded sedations.**

In the MENU under DOCUMENTS it is possible to:

- view a recorded sedation by choosing from the list displayed in chronological order
- export the recorded sedations to a USB stick, after inserting the stick into the dedicated socket, press the appropriate button, select whether to export the last 99 sedations or all recorded sedations and press export.

Dedicated software is required to view the exported data.

**Exported data is not deleted from the device.**

**The file name is sorted by sequential numbers e.g. MF000001**

**NOTE:**

**Tecno-gaz declines all responsibility in the event of data loss. Ensuring the security and protection of patient data/sedation protocols is the responsibility of the operator in charge of the device.**

**IT IS RECOMMENDED TO EXPORT THE DATA TO A USB STICK AND TO STORE IT IN A SECURE DATABASE.**

## 7 ALARMS

### 7.1 WARNINGS

An alarm in the system is a recall caused by system problems that are sent to medical personnel by means of a VISUAL AND AUDIO MESSAGE, when the delivery pressure in the equipment is insufficient or a device failure occurs that may compromise the safety of the device.

	THE DEVICE DOES NOT MONITOR THE PATIENT'S PHYSIOLOGICAL PARAMETERS. DURING SEDATION MONITOR THE PATIENT USING A PULSE OXIMETER WITH CUSTOMISED ALARM THRESHOLDS.
	AT START-UP, TO CHECK THE OPERATION OF THE ACOUSTIC SIGNAL, THE DEVICE EMITS 3 BEEPS, IN CASE OF ANOMALY DO NOT USE THE DEVICE AND CONTACT THE TECHNICAL SUPPORT.

The device warns the operator by means of WARNINGS if there are any anomalies, THESE ANOMALIES DO NOT AFFECT THE SAFETY OF THE DEVICE.

The WARNING indication appears on the screen, accompanied by an acoustic alarm that repeats with 3 beeps and a pause.

If the device is connected to the mains, the BARLED (see par 5.2) changes to Yellow.

To continue using the device it is necessary to press "Reset", the WARNING screen and the BUZZER SIGNAL are removed, a BLINKING ICON remains on the DISPLAY and the BARLED is yellow.

**If the cause of the fault is removed, the icon on the display disappears.**

**Below is a list of warnings and possible remedial actions:**

Warning	Behaviour	Possible cause	Possible solution
A1*	The device does not deliver gas, the O <sub>2</sub> pressure measured is 0	Closed or empty O <sub>2</sub> cylinder	Open the O <sub>2</sub> cylinder Replace the O <sub>2</sub> cylinder
A2*	The device does not deliver N <sub>2</sub> O, N <sub>2</sub> O pressure measured equal to 0	Closed or empty N <sub>2</sub> O cylinder	Open the N <sub>2</sub> O cylinder Replace the N <sub>2</sub> O cylinder
A3*	Device only delivers O <sub>2</sub> , low O <sub>2</sub> pressure measured (< 2 bar)	Possible problem with the O <sub>2</sub> reducer	Open the O <sub>2</sub> cylinder Replace the O <sub>2</sub> cylinder
A5*	The device does not deliver gas, high O <sub>2</sub> pressure measured (> 6 bar)		Check the O <sub>2</sub> pressure reducer **
A4*	Device only delivers O <sub>2</sub> , low N <sub>2</sub> O pressure measured (< 2 bar)	Possible problem with N <sub>2</sub> O reducer	Open the N <sub>2</sub> O cylinder Replace the N <sub>2</sub> O cylinder
A6*	Device only delivers O <sub>2</sub> , high N <sub>2</sub> O pressure measured (> 6 bar)		Check the N <sub>2</sub> O pressure reducer **
A7***	Battery < 20% The device signals that the battery is running low.	-	Connect the device to the mains
A8	Battery ≤ 0% The device interrupts the cycle Low battery	-	Connect the device to the mains or switch off the device to avoid damaging the battery

\* This warning appears only during sedation

\*\* If the pressure relief valve on the reducer vents, the reducer must be replaced.

\*\*\* Only if a sedation is being performed

## 7.2 - ERRORS

The device warns the operator by means of ERRORS if there are any anomalies, THESE ANOMALIES AFFECT THE SAFETY OF THE DEVICE.

The ERROR indication appears on the screen, accompanied by an acoustic alarm that repeats with 3 beeps, pause and 2 beeps.

If the device is connected to the mains, BARLED (see par 5.1) changes to Red.

Error	Behaviour	Possible cause
E1	O <sub>2</sub> loss	O <sub>2</sub> solenoid valve leakage
E2	N <sub>2</sub> O loss	N <sub>2</sub> O solenoid valve leakage
E4*	Inequality between O <sub>2</sub> flow delivered and flow set	Incorrect O <sub>2</sub> knob alignment. Flowmeter problems.
E5*	Inequality in N <sub>2</sub> O flow delivered and flow set	N <sub>2</sub> O knob alignment not correct. Flowmeter problems.
E7*	% nitrous oxide greater than 70	Possible calibration error
E8	O <sub>2</sub> pressure sensor problems	Pressure sensor breakage Damaged wiring
E9	N <sub>2</sub> O pressure sensor problems	Pressure sensor breakage Damaged wiring
E10	O <sub>2</sub> flowmeter problems	Flowmeter breakage Damaged wiring
E11	N <sub>2</sub> O flowmeter problems	Flowmeter breakage Damaged wiring

\* This warning appears only during sedation

If during sedation the device detects a flow rate higher than 70% (maximum 1% within the tolerances related to the mechanics of the device) the device warns the operator with the activation of the acoustic signal and the indication of the % turns red and flashes for 10 sec, if the percentage is reduced it is possible to terminate sedation, if the percentage is not reduced the device stops the delivery and the ERROR sign appears.

IT IS POSSIBLE TO RESET THE ERROR BY ENTERING THE PASSWORD: 1234

In any case, it is necessary to contact technical support.

**IN THE EVENT OF AN ERROR, DO NOT USE THE DEVICE  
SWITCH OFF THE DEVICE AND CONTACT THE TECHNICAL SUPPORT**

## 7.3 WI-FI COMMUNICATION

The device alerts the operator, if the device is connected, on the status of the WI-FI connection and data sending.

A number from 1 to 8, next to the WI-FI icon, identifies the type of problem.

**If there are any problems, CHECK THAT THE MAINS NAME AND PASSWORD ARE ENTERED CORRECTLY.**

If the problem persists, contact technical support.

## 8 MAINTENANCE

### BEFORE EACH USE, CHECK

THE DELIVERY: Start the delivery and check that there are no anomalies.

N2O BLOCK: Set a flow rate of 8 l/min with 50% N2O concentration, close the oxygen cylinder (simulating a finished cylinder situation) and make sure that after a few seconds, the N2O supply is completely blocked; if only N2O is delivered, do not use the device and contact technical support.

FLUSH BUTTON: Set up a supply with 50% concentration, check that the nitrous oxide supply is interrupted and the oxygen flow is at maximum.

### EVERY WEEK

#### CHECK FOR LEAKS:

- Open the valves of the Oxygen and Nitrous Oxide cylinders and note the values indicated by the pressure gauges. Close cylinder valves. Wait ~ 5 minutes and measure the new values indicated by the gauges. If the values coincide, the device does not leak.

- Disconnect the delivery tube from the mask, deliver approx. 5 l/min O<sub>2</sub>, close the tube with one hand, check that the anaesthesia bag is inflated and that there are no leaks (e.g. tube cut, punctured balloon).

Interrupt the flow.

### 8.1 ROUTINE MAINTENANCE - CLEANING

The DEVICE must be cleaned using clean gloves in a clean environment.

Use a disinfectant that is compatible with plastics (avoid cream products).

Do not spray the disinfectant directly onto the component to be cleaned, but onto a cloth.

Use disinfectants in the dilution indicated by the manufacturer.

The cleaning processes for individual components are described in the enclosed cleaning manual (e.g. circuits, masks, reducers, etc.)

### 8.2 ROUTINE MAINTENANCE - CYLINDER REPLACEMENT

Cylinders must be replaced with extreme caution (see 5.5 CYLINDERS).

Procedure:

1. Slowly close the valves (turn CLOCKWISE).
2. If the system is under pressure, relieve the pressure to make it easier to remove the reducers.
3. Store the reducers in a sheltered, clean place.
4. Install the reducers on the filled cylinders.
5. OPEN THE VALVES SLOWLY, A RAPID OPENING OF THE CYLINDER VALVE CAN CAUSE A SUDDEN RISE IN TEMPERATURE INSIDE THE REDUCER, WHICH CAN LEAD TO DAMAGE AND BREAKAGE.

**Refilling, inspection and periodic testing of cylinders must be carried out by specialised companies in accordance with the directives in force in the country of use.**

### 8.3 EXTRAORDINARY MAINTENANCE

THE DEVICE MUST UNDERGO A FUNCTIONAL AND ELECTRICAL SAFETY CHECK (ELECTROMEDICAL EN 62353).

THE RECOMMENDED FREQUENCY OF SAFETY INSPECTIONS IS ONCE A YEAR FOR EQUIPMENT USED IN SURGICAL AND SIMILAR PREMISES AND ONCE EVERY TWO YEARS FOR EQUIPMENT USED IN ALL OTHER TYPES OF PREMISES SUCH AS DENTAL PRACTICES.

THE ELECTRICAL SAFETY CHECK MUST ALWAYS BE CARRIED OUT AFTER REPAIR WORK ON THE EQUIPMENT.

CONTACT TECHNICAL SUPPORT TO BOOK AN INTERVENTION.

### 8.4 PROBLEMS - FAILURE TO SEDATE THE PATIENT

Possible causes	Solution
The mask does not adhere to the patient's face.	Use a mask with a size suitable for the patient.
The correct % of N <sub>2</sub> O was not determined.	BASELINE is characteristic of the patient, it is important to assess the characteristic parameters that indicate that the patient is sedated.
The one-way valves on the mask are damaged or the membrane is missing.	Replace the mask.
The patient does not sedate.	Instruct the patient to breathe through the nose.
Active suction is excessive and takes up the gas before it is taken up by the patient.	Check, the behaviour of the anaesthesia bag, it must remain open and follow the patient's breathing If necessary, reduce the suction flow.
The circuit and/or the anaesthesia bag are damaged.	Small cuts or cracks may allow gas to escape and not reach the patient. Replace the damaged component.
Leakage of the emergency suction valve.	Set an oxygen flow. Close the delivery hose connection (see paragraph 6.1 position A) with one hand and check that the balloon inflates. Remove your hand and stop the flow. If any leaks are noticed, contact the Tecno-gaz service centre.
The anaesthesia bag inflates until it detaches from the hose connection.	The gas supply pipe to the patient was reversed with the exhaled gas pipe

## 9 SUPPORT

Technical support must be requested from the depot that invoiced the device, or directly to:

TECNO-GAZ S.p.A. Strada Cavalli, 4 – 43038 - Sala Baganza Parma ITALY. www.tecnogaz.com	PHONE +39 0521 8380 FAX +39 0521 833391 support@tecnogaz.com
--	--

The technical support centre will evaluate the return to the premises or the intervention of a technician and, having inspected the machine, will draw up an estimate of costs, which will be forwarded to the distributor customer who will transmit it to the end customer for acknowledgement and signature.

Upon receipt of the signed quotation for acceptance, the device will be put into production and will be shipped within the time specified on the quotation form.

For shipping the device:

- Use the original packaging; if the original packaging is no longer available, use suitable packaging. The goods travel at the sender's risk.
- Ship the complete device.

The reducers must not be contaminated (detergents, alcohol, oils ...). If the technical support department does not consider the packaging to be suitable, the reducers will be scrapped.

- Indicate in writing, and include in the packaging, a document indicating precisely the anomaly detected or the service intend to be used.
- Always ship carriage paid, otherwise the transport costs incurred will be charged.

All non-original packaging received will be disposed of.

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