

ORION TECH

Surgical light




Orion Tech is a surgical light designed to meet the needs and overcome the limitations of dental surgery:

- A unique field concept specially created to illuminate the mouth in any position, including when angled at 45° both horizontally and vertically.
- Two colour temperatures to distinguish between soft tissue, hard tissue and blood.
- Wide range of light intensity variation taking into account the reflection of white surfaces.
- Unique, patented indirect light concept ensures constant intensity and reduced heat emission.

Stable and uniform light

Patented indirect lighting technology eliminates shadows, offering unprecedented illumination.



Orion Tech Surgical light



Quality that lasts over time

Light source

Protected light source
resistant to air, water and dust

Dome



Functionality and design

Minimal lines and cutting-edge shapes

Hygiene and ease of cleaning




Illuminated fields

Elliptical field 20x15 cm
For dental treatments


Circular field 19 cm and 27 cm
For surgical procedures




SINGLE FORK VERSION



DOUBLE FORK VERSION



LC003LRE
ON TROLLEY
VERSION



LC001LRE
CEILING DOUBLE
FORK VERSION



LC008LRE
CEILING SINGLE
FORK VERSION



LC005LRE
CEILING VERSION
WITH DOUBLE LAMP



LC006LRE
CEILING WITH
CAMERA SUPPORT